What's new?

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Korean Society of Gastrointestinal Endoscopy (KSGE) Guidelines for Endoscope Reprocessing

Byung-Wook Kim¹, Byoung Kwan Son², Won Hee Kim³, Dae-Sung Myung⁴, Young-Seok Cho⁵, and Byung Ik Jang⁶;
The Disinfection Management and Conscious Sedation Committee, Korean Society of Gastrointestinal Endoscopy

¹*Division of Gastroenterology, Department of Internal Medicine, Incheon St. Mary's Hospital, College of Medicine, The Catholic University of Korea, Incheon;* ²*Department of Internal Medicine, Eulji Hospital, Eulji University School of Medicine, Seoul;* ³*Digestive Disease Center, CHA Bundang Medical Center, CHA University, Seongnam;* ⁴*Departments of Internal Medicine, Chonnam National University Medical School, Gwangju;* ⁵*Department of Internal Medicine, College of Medicine, The Catholic University of Korea, Seoul;* ⁶*Department of Internal Medicine, Yeungnam University College of Medicine, Daegu, Korea*

Korean Society of Gastrointestinal Endoscopy (KSGE) issued its guidelines for endoscope reprocessing for the first time in 1995, and the version of the guidelines has been updated in August 2009, August 2012 and March 2015, respectively. Guidelines for endoscope reprocessing should be revised continuously because new disinfectants and devices are introduced. The current official version of the KSGE guidelines for endoscope reprocessing is explained herein to assist the reader in understanding of KSGE requirements for cleaning and disinfection of endoscopes. (Korean J Med 2017;92:239-244)

Keywords: Endoscopy; Guidelines; Endoscope reprocessing

Correspondence to Young-Seok Cho, M.D.
Department of Internal Medicine, Seoul St. Mary's Hospital, College of Medicine, The Catholic University of Korea, 222 Banpo-daero, Seocho-gu, Seoul 06591, Korea
Tel: +82-2-2258-6021, Fax: +82-2-2258-2038, E-mail: yscho@catholic.ac.kr

*This manuscript has been published in Clinical Endoscopy (2017;50:143-147) as English, but also published in The Korean Journal of Medicine as Korean to broaden distribution of its contents.

Copyright © 2017 The Korean Association of Internal Medicine
This is an Open Access article distributed under the terms of the Creative Commons Attribution Non-Commercial License (<http://creativecommons.org/licenses/by-nc/3.0/>) which permits unrestricted noncommercial use, distribution, and reproduction in any medium, provided the original work is properly cited.

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1. Precleaning

- (1) Upon completion of endoscopy, immediately clean the exterior of the endoscope beside the bed using gauze soaked in enzymatic detergent or sterilized distilled water.
- (2) Place the distal end of the endoscope into detergent solution, aspirate detergent through the air/water and biopsy channels and then blow air through. These steps should be repeated to remove residual contamination
- (3) After disconnecting the endoscope from the power source, transfer to a separately located reprocessing area in a designated container. If the distance to the reprocessing area is far, the container should be enclosed to prevent a possible infection.

2. Cleaning

- (1) Remove all the detachable parts from the endoscope
- (2) Detect if there is a leak.
- (3) Clean the endoscope with detergent solution and thoroughly brush all channels and removable parts.
- (4) Clean the parts that cannot adequately brushed using an ultrasonic cleaner while they are immersed in detergent solution.
- (5) Thoroughly rinse the endoscope and detachable parts with clean water to wash out residual detergent (even if in case of using an automatic reprocessor, cleaning needs to be done manually.)

3. Disinfection

- (1) In case of manual disinfection, use a basin that is large enough and corrosion-resistant. Completely immerse the endoscope and all removable parts completely in a basin of high-level disinfectants; ensure that all channels are completely filled with disinfectant.
- (2) When using any disinfectant solution, follow manufacturer's instructions for disinfection conditions and time.
- (3) If an automatic reprocessor is used, follow all instructions provided in User's Guide.

4. Rinse

- (1) Thoroughly rinse the endoscope and flush all channels using clean drinkable water.

5. Drying

- (1) Flush each channel with compressed air and 70 to 90% ethyl or isopropyl alcohol to facilitate drying.

6. Storage

- (1) Store endoscopes in a designated, well-ventilated area. By hanging endoscopes in a vertical position, ensure that their distal end do not touch the floor.

7. Endoscopic accessories

- (1) All reusable endoscopic accessories should be disinfected.
- (2) Immerse dismantled endoscopic accessories in detergent solution and use brush and sponge to clean the internal surface.
- (3) Clean accessories using an ultrasonic cleaner in detergent solution.
- (4) After rinsing accessories with clean water and dry with clean cloth and compressed air.
- (5) Before storage, sterilize or disinfect accessories, according to the type of accessories and the corresponding manufacturer's instructions. Reusable accessories such as suction, air water valve and rubber biopsy valve caps should be high-level disinfected before reuse.
- (6) Disposable needles or forceps should not be used more than once.

8. Water bottle and connecting tube

- (1) The water bottle and its connecting tube should be sterilized once a day and contain sterile water.

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REFERENCES

1. Disinfection Management Committee, Korean Society of Gastrointestinal Endoscopy. Guidelines for disinfection of gastrointestinal endoscopes. Seoul: Medbooks, 2012:23-33.
2. Ministry of Health and Welfare. Guidelines for medical instruments and supplies. Notification No. 2010-61, Ministry of Health and Welfare. Sejong: Ministry of Health and Welfare, 2010.
3. Centers for Disease Control and Prevention. Sterilization and disinfection guidelines in the medical institution. Cheongju: Centers for Disease Control and Prevention, 2014.
4. Oh HJ, Kim JS. Clinical practice guidelines for endoscope reprocessing. Clin Endosc 2015;48:364-368.
5. Hookey L, Armstrong D, Enns R, Matlow A, Singh H, Love J. Summary of guidelines for infection prevention and control for flexible gastrointestinal endoscopy. Can J

- Gastroenterol 2013;27:347-350.
6. ASGE Quality Assurance In Endoscopy Committee, Petersen BT, Chennat J, et al. Multisociety guideline on reprocessing flexible gastrointestinal endoscopes: 2011. Gastrointest Endosc 2011;73:1075-1084.
 7. Society of Gastroenterology Nurses and Associates, Inc. SGNA standards: standards of infection control in reprocessing of flexible gastrointestinal endoscopes. Gastroenterol Nurs 2010;33:70-80.
 8. Yoon IJ. Tips for safe endoscopic reprocessing. 49th Seminar of Korean Society of Gastrointestinal Endoscopy Proceedings. 2013 Aug 25; Goyang. 2013:324-327.
 9. Son BK. Disinfectants for endoscope reprocessing. 50th Seminar of Korean Society of Gastrointestinal Endoscopy Proceedings. 2014 Mar 30; Goyang. 2014:327-329.
 10. Oh KJ. Reprocessing and maintaining of endoscope accessories. 50th Seminar of Korean Society of Gastrointestinal Endoscopy Proceedings. 2014 Mar 30; Goyang. 2014:330-334.
 11. Cho YS. Importance of endoscope reprocessing and new guideline for reprocessing. 52nd Seminar of Korean Society of Gastrointestinal Endoscopy Proceedings. 2015 Mar 29; Goyang. 2015:338-342.
 12. Lim CH, Choi MG, Kim WC, et al. Performance and cost of disposable biopsy forceps in upper gastrointestinal endoscopy: comparison with reusable biopsy forceps. Clin Endosc 2012;45:62-66.